

Section 1: Identification of the Substance/Mixture and of the Company/Undertaking

Product identifier:

Identification as on the label/Trade name: Morris Reagent Part B

Relevant identification uses of the substance and uses advised against:

Identified uses: TN Scientific's Morris Reagent is designed for the qualitative detection of benzodiazepines in substances, intended for laboratory and professional use to aid in substance identification and harm reduction.

Uses advised against: No other uses are advised.

Details of the supplier of the Safety Data Sheet:

TN Scientific LLC
6700 Baum Drive Suite 18
Knoxville, TN
37931 USA
+1-865-243-6011

Emergency telephone numbers:

24-hour Emergency Contact:
CHEMTREC 24-hour: +1-800-424-9300

Section 2: Hazards Identification

Classification of the substances or mixture:

The mixture is classified according to: Regulation EC 1272/2008 [EU-GHS/CLP]

Hazard classes/Hazard categories:

Skin irritant (Category 2)
Eye irritant (Category 2)

Label elements:

Hazard pictogram:



Signal word:Warning.

Hazard statements:

H315 Causes skin irritation.
H319 Causes serious eye irritation.

Precautionary statements:

Prevention

P264 + P265 Wash hands thoroughly after handling. Do not touch eyes.
P280 Wear protective gloves/eye protection.

Response

P302 + P352 IF ON SKIN: Wash with plenty of water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P332 + P317 If skin irritation occurs: Get medical help.

P337 + P317 If eye irritation persists: Get medical help.

P362 + P364 Take off contaminated clothing and wash it before reuse.

Section 4: First-Aid Measures

Description of first aid measures:

If inhaled: If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.

In case of skin contact: Wash off with soap and plenty of water. Take victim immediately to hospital. Consult a physician.

In case of eye contact: Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.

If swallowed: Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.

Most important symptoms and effects, both acute and delayed:

Causes skin and eye irritation.

Indication of any immediate medical attention and special treatment needed: Treat symptomatically.

Section 5: Fire-Fighting Measures

Extinguisher media:

Suitable extinguishing media: Dry chemical, CO₂, water spray or alcohol-resistant foam.

Unsuitable extinguishing media: None known.

Special hazards arising from the substance or mixture: None known.

Advice for fire-fighters: As in any fire, wear self-contained breathing apparatus and full protective gear.

Section 6: Accidental Release Measures

Personal precautions, protective equipment and emergency procedures:

Avoid breathing vapors. Ensure adequate ventilation. Evacuate personnel to safe areas.

Environmental precautions: Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

Methods for containment and cleaning up:

Wipe up spills with absorbent material. If necessary clean with disinfectant afterwards and dispose of in accordance with the local regulations. Rinse affected area with water afterwards.

Section 7: Handling and Storage

Precautions for safe handling:

Handle in accordance with good industrial hygiene and safety practice. Avoid contact with skin and eyes.

Conditions for safe storage, including incompatibilities:

Store in cool, dry place. Keep container tightly closed.

Section 8: Exposure Controls and Personal Protection

Control parameters:

Occupational exposure limits: Does not contain any hazardous materials with Occupational Exposure Limits.

Exposure controls:

Appropriate engineering controls: Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the product.

Individual protection measures, such as personal protective equipment:

Eye/face protection: Safety glasses.

Skin protection: Handle with gloves. Gloves must be inspected prior to use. Use proper glove removal technique (without touching glove's outer surface) to avoid skin contact with this product. Dispose of contaminated gloves after use in accordance with applicable laws and good laboratory practices.

Body protection: Protective clothing as appropriate.

Respiratory protection: No respiratory protection is required.

Section 9: Physical and Chemical Properties

Information on basic physical and chemical properties:

Appearance (form): Liquid.

Color: Colorless.

Odor: Odorless.

Odor threshold: No data available.

pH (concentration): No data available.

Melting point/range (°C): 0 °C

Boiling point/range (°C): 100 °C

Flash point (°C): No data available.

Evaporation rate: No data available. **Flammability (solid, gas):** No data available. **Upper/lower flammability/explosive limits:** No data available. **Vapor pressure:** No data available. **Vapor density:** No data available. **Relative density:** 1.0 **Water solubility (g/L):** Soluble. **n-Octanol/Water partition coefficient:** No data available. **Auto-ignition temperature:** No data available. **Viscosity, dynamic (mPa.s):** No data available.

Section 10: Stability and Reactivity

Reactivity: No specific test data related to reactivity available for this product or its ingredients.
Chemical stability: Stable under recommended storage conditions.
Possibility of hazardous reactions: No hazardous reactions are known.
Conditions to avoid: None known.
Incompatible materials: None known.
Hazardous decomposition products: Under normal conditions of storage and use, hazardous decomposition products should not be produced.

Section 11: Toxicological Information

Information on toxicological effects:

Acute toxicity: No data available.
Skin corrosion/irritation: Causes skin irritation.
Serious eye damage/irritation: Causes serious eye irritation.
Respiratory or skin sensitization: No data available.
Germ cell mutagenicity: No data available.
Carcinogenicity: No data available.
Reproductive toxicity: No data available.
STOT-single exposure: No data available.
STOT-repeated exposure: No data available.
Aspiration hazard: No data available.

Section 12: Ecological Information

Toxicity: No data available.
Persistence and degradability: No data available.
Bioaccumulative potential: No data available.
Mobility in soil: No data available.
Results of PBT& vPvB assessment: No data available.

Section 13: Disposal Considerations

Waste treatment methods: Dispose of in accordance with local regulations.
Product/packaging disposal: Dispose of as unused product.

Section 14: Transport Information

UN number: -

UN proper shipping name: Not regulated as dangerous goods.

Transport hazard class: -

Packing group: -

Environmental hazards: -

Special precautions for user: Not dangerous goods in the meaning of ADR/RID, ADN, IMDG-Code, ICAO/IATA-DGR.

Transport in bulk according to Annex II of Marpol and the IBC Code: Not applicable.

Section 15: Regulatory Information

Safety, health and environmental regulations/legislation for the mixture: EU Regulations: Regulation EC 1272/2008 [EU-GHS/CLP]

US Regulations:

SARA Title III Section 302/304 Extremely Hazardous Substance: No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA Title III Section 313 Supplier Information: This material does not contain any chemical components with known CAS numbers that exceed the threshold (De Minimis) reporting levels established by SARA Title III, Section 313.

CERCLA Section 102(a) Hazardous Substance: Sodium hydroxide, CAS No. 1310-73-2: 1,000 lbs.

California Proposition 65: This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

Chemical Safety Assessment carried out: No.

Section 16: Other Information

Indication of changes: GHS aligned.

Relevant classification and H statements (number and full text):

H314 Causes severe skin burns and eye damage.

H315 Cause skin irritation.

H319 Causes serious eye irritation.

NFPA rating:



Health: 2

Fire: 0

Reactivity: 0

Training instructions: Use as instructed.

Further information: This information is based upon the present state of our knowledge. This SDS has been compiled and is solely intended for this product.

Notice to readers: Employers should use this information only as a supplement to other information gathered by them and should make independent judgment of suitability of this information to ensure proper use and protect the health and safety of employees.



This information is furnished without warranty, and any use of the product not in conformance with this Safety Data Sheet, or in combination with any other product or process, is the responsibility of the user.